


MYF5 Antibody (R34985)

Catalog No.	Formulation	Size
R34985-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4617
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MYF5 antibody is available for research use only.



Western blot testing of muscle lysate with MYF5 antibody at 0.5ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight ~28kDa.

Description

Additional name(s) for this target protein: Myogenic factor 5

Application Notes

Optimal dilution of the MYF5 antibody should be determined by the researcher.

Immunogen

Amino acids PECNSPVWSRKSST were used as the immunogen for this MYF5 antibody.

Storage

Aliquot and store the MYF5 antibody at -20oC.