


MYH1 Antibody [MUTYH] (R34210)

Catalog No.	Formulation	Size
R34210-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4595
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This MYH1 antibody is available for research use only.



Western blot testing of human bone marrow lysate with MYH1 antibody at 0.1ug/ml.
Predicted molecular weight: ~60kDa.

Description

Additional name(s) for this target protein: MutY homolog 1

Application Notes

Optimal dilution of the MYH1 antibody should be determined by the researcher.

Immunogen

Amino acids HISTDAHSLNSAAQ were used as the immunogen for this MYH1 antibody.

Storage

Aliquot and store the MYH1 antibody at -20oC.