


MYH9 Antibody (R33596)


Catalog No.	Formulation	Size
R33596-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4627
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MYH9 antibody is available for research use only.



Western blot testing of human peripheral blood lymphocyte lysate with MYH9 antibody at 0.3ug/ml. Predicted molecular weight: ~227kDa.



Western blot testing of NIH3T3 lysate with MYH9 antibody at 1ug/ml. Predicted molecular weight: ~227kDa.

Description

Additional name(s) for this target protein: Myosin, heavy chain 9

Application Notes

Optimal dilution of the MYH9 antibody should be determined by the researcher.

Immunogen

Amino acids DQINTDLNLERSH were used as the immunogen for this MYH9 antibody.

Storage

Aliquot and store the MYH9 antibody at -20oC.