


MYLK Antibody / Myosin light chain kinase (R34387)

Catalog No.	Formulation	Size
R34387-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4638
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This MYLK antibody is available for research use only.



Description

Additional name(s) for this target protein: Myosin, light chain kinase

Application Notes

Optimal dilution of the MYLK antibody should be determined by the researcher.

Immunogen

Amino acids KDTKNMEAKKLSKD were used as the immunogen for this MYLK antibody. **Storage** Aliquot and store the MYLK antibody at -20oC.