

MYO9B Antibody (R34539)

Catalog No.	Formulation	Size
R34539-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4650
Applications	Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This MYO9B antibody is available for research use only.



Western blot testing of human liver lysate with MYO9B antibody at 2ug/ml. Predicted molecular weight: ~243kDa, observed here at ~200kDa.

Description

Additional name(s) for this target protein: Myosin IXB

Application Notes

Optimal dilution of the MYO9B antibody should be determined by the researcher.

Immunogen

Amino acids KNDFEPVKQSKAQ were used as the immunogen for this MYO9B antibody.

Storage

Aliquot and store the MYO9B antibody at -20oC.