

Myogenin Antibody (R36140)

Catalog No.	Formulation	Size
R36140-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P15173
Gene ID	4656
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Myogenin antibody is available for research use only.



Western blot testing of mouse skeletal muscle lysate with Myogenin antibody at 2ug/ml.
Molecular weight: routinely observed at ~34 kDa.

Description

Additional name(s) for this target protein: MYOG, MYF4, myogenic factor 4; MYOG

Application Notes

Optimal dilution of the Myogenin antibody should be determined by the researcher.

Immunogen

Amino acids TDAHNLHSLTSIVD were used as the immunogen for this Myogenin antibody.

Storage

Aliquot and store the Myogenin antibody at -20°C.