


NAC1 Antibody / BTBD14B (R34663)

Catalog No.	Formulation	Size
R34663-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	112939
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This NAC1 antibody is available for research use only.



Western blot testing of human bone marrow lysate with NAC1 antibody at 0.5ug/ml.
Predicted molecular weight: ~57 kDa, commonly observed at 57-70 kDa.

Description

Additional name(s) for this target protein: BTB (POZ) domain containing 14B; BTBD14B

Application Notes

Optimal dilution of the NAC1 antibody should be determined by the researcher.

Immunogen

Amino acids KTEQQESDSVQC were used as the immunogen for this NAC1 antibody.

Storage

Aliquot and store the NAC1 antibody at -20oC.