


## NALP3 Antibody / CIAS1 (R33843)

Catalog No.	Formulation	Size
R33843-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114548
Applications	Flow Cytometry : 10ug/10 <sup>6</sup> cells in 0.2ml ELISA (peptide) LOD : 1:32000
Limitations	This NALP3 antibody is available for research use only.



FACS testing of fixed and permeabilized human U937 cells with NALP3 antibody (blue) at 10ug/10<sup>6</sup> cells and naive goat Ig (black).

## Description

Additional name(s) for this target protein: NACHT LRR and PYD containing protein 3; NLRP3; CIAS1

## Application Notes

Optimal dilution of the NALP3 antibody should be determined by the researcher.

1. Reported to detect NALP3 by FACS in human monocytes

## **Immunogen**

Amino acids EKPELTVVFEPWS were used as the immunogen for this NALP3 antibody.

## **Storage**

Aliquot and store the NALP3 antibody at -20°C.