


NANOS Antibody / NANOS1 (R35839)

Catalog No.	Formulation	Size
R35839-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	340719
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This NANOS antibody is available for research use only.



Western blot testing of human testis lysate with NANOS antibody at 0.3ug/ml. Predicted molecular weight: ~30kDa but routinely observed at 35~40 kDa (see ref. 1).

Description

Additional name(s) for this target protein: NOS1, NANOS1

Application Notes

Optimal dilution of the NANOS antibody should be determined by the researcher.

Immunogen

Amino acids SARDGPPGKKLR were used as the immunogen for this NANOS antibody.

Storage

Aliquot and store the NANOS antibody at -20oC.

References (1)