

NAT1 Antibody (R34827)

Catalog No.	Formulation	Size
R34827-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NAT1 antibody is available for research use only.



Western blot testing of human erythrocytes lysate with NAT1 antibody at 1ug/ml.
Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: N-acetyltransferase 1

Application Notes

Optimal dilution of the NAT1 antibody should be determined by the researcher.

Immunogen

Amino acids NISLQRKLVPKHGDR were used as the immunogen for this NAT1 antibody.

Storage

Aliquot and store the NAT1 antibody at -20oC.