

## NCAPH Antibody (R36151)

Catalog No.	Formulation	Size
R36151-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23397
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NCAPH antibody is available for research use only.



Western blot testing of K562 cell lysate with NCAPH antibody at 0.01ug/ml. Predicted molecular weight: ~83kDa, routinely observed at 90~97kDa.

## Description

Additional name(s) for this target protein: Non-SMC condensin I complex, subunit H; HCAP-H

## Application Notes

Optimal dilution of the NCAPH antibody should be determined by the researcher.

## Immunogen

Amino acids TEDLSDVLVRQGD were used as the immunogen for this NCAPH antibody.

## Storage

Aliquot and store the NCAPH antibody at -20oC.

## References (1)