


NDUFA7 Antibody (R35825)

Catalog No.	Formulation	Size
R35825-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4701
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NDUFA7 antibody is available for research use only.



Western blot testing of human heart lysate with NDUFA7 antibody at 0.1 ug/ml. Predicted molecular weight: ~13kDa.

Description

Additional name(s) for this target protein: NADH dehydrogenase (ubiquinone) 1 alpha subcomplex 7

Application Notes

Optimal dilution of the NDUFA7 antibody should be determined by the researcher.

Immunogen

Amino acids RYQEISKRTQPP were used as the immunogen for this NDUFA7 antibody.

Storage

Aliquot and store the NDUFA7 antibody at -20oC.