

NDUFS1 Antibody (R36160)

Catalog No.	Formulation	Size
R36160-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rabbit, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4719
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This NDUFS1 antibody is available for research use only.



Western blot testing of human heart lysate with NDUFS1 antibody at 0.3ug/ml. Predicted molecular weight: ~80kDa.

Description

Additional name(s) for this target protein: NADH dehydrogenase (ubiquinone) Fe-S protein 1

Application Notes

Optimal dilution of the NDUFS1 antibody should be determined by the researcher.

Immunogen

Amino acids TEKSATYVNTEGR were used as the immunogen for this NDUFS1 antibody.

Storage

Aliquot and store the NDUFS1 antibody at -20oC.