


NDUFS3 Antibody (R33662)

Catalog No.	Formulation	Size
R33662-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4722
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This NDUFS3 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with NDUFS3 antibody at 0.3ug/ml.
Predicted molecular weight: ~30 kDa.

Description

Additional name(s) for this target protein: NADH dehydrogenase (ubiquinone) Fe-S protein 3

Application Notes

Optimal dilution of the NDUFS3 antibody should be determined by the researcher.

Immunogen

Amino acids DTRPTVRPRNDVAHK were used as the immunogen for this NDUFS3 antibody.

Storage

Aliquot and store the NDUFS3 antibody at -20oC.