


NDUFS3 Antibody (R33662)

Catalog No.	Formulation	Size
R33662-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4722
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This NDUFS3 antibody is available for research use only.



Description

Additional name(s) for this target protein: NADH dehydrogenase (ubiquinone) Fe-S protein 3

Application Notes

Optimal dilution of the NDUFS3 antibody should be determined by the researcher.

Immunogen

Amino acids DTRPTVRPRNDVAHK were used as the immunogen for this NDUFS3 antibody. **Storage** Aliquot and store the NDUFS3 antibody at -20oC.