


NEDD8 Antibody (R32590)


Catalog No.	Formulation	Size
R32590	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
UniProt	Q15843
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml Flow Cytometry : 1-3ug/10 ⁶ cells Immunofluorescence (FFPE) : 2-4ug/ml
Limitations	This NEDD8 antibody is available for research use only.




Immunofluorescent staining of FFPE human A431 cells with NEDD8 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



IHC testing of FFPE human lung cancer tissue with NEDD8 antibody at 1ug/ml. HIER: steam section in pH6 citrate buffer for 20 min.



IHC testing of FFPE human lung cancer tissue with NEDD8 antibody at 1ug/ml. HIER: steam section in pH6 citrate buffer for 20 min.



Flow cytometry testing of human A431 cells with NEDD8 antibody at 1ug/10⁶ cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= NEDD8 antibody.



Western blot testing of 1) rat testis, 2) mouse thymus, 3) mouse brain, 4) human HeLa and 5) human MCF7 lysate with NEDD8 antibody at 0.5ug/ml. Predicted molecular weight ~9 kDa.

Description

NEDD8 is a protein that in humans is encoded by the NEDD8 gene. Human NEDD8 shares 60% amino acid sequence identity to ubiquitin. The most substrates of NEDD8 modification are the Cullin subunits of Cullin-based E3 ubiquitin ligases, which are active only when neddylated. Their NEDDylation is critical for the recruitment of E2 to the ligase complex, thus facilitating ubiquitin conjugation. NEDD8 modification has therefore been implicated in cell cycle progression and cytoskeletal regulation.

Application Notes

Optimal dilution of the NEDD8 antibody should be determined by the researcher.

Immunogen

Amino acids 20-60 (TDKVERIKERVEEKEGIPPQQRLIYSGKQMNDEKTAADYK-human) were used as the immunogen for the NEDD8 antibody.

Storage

After reconstitution, the NEDD8 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at

-20oC. Avoid repeated freezing and thawing.