

NEDD9 Antibody / hEF1 / CAS-L (R34845)

Catalog No.	Formulation	Size
R34845-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4739
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This NEDD9 antibody is available for research use only.



Western blot testing of human kidney lysate with NEDD9 antibody at 0.1ug/ml. Predicted molecular weight: ~93kDa.

Description

Additional name(s) for this target protein: Neural precursor cell expressed, developmentally down-regulated 9, CASL, HEF1; NEDD9; Enhancer of filamentation 1

Application Notes

Optimal dilution of the NEDD9 antibody should be determined by the researcher.

1. This NEDD9 antibody is specific for isoform 1.

Immunogen

Amino acids NKPQNKDDLDR were used as the immunogen for this NEDD9 antibody.

Storage

Aliquot and store the NEDD9 antibody at -20oC.