


NEK7 Antibody (R33252)


Catalog No.	Formulation	Size
R33252-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	140609
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NEK7 antibody is available for research use only.



Western blot testing of human lung lysate with NEK7 antibody at 0.1ug/ml. Predicted molecular weight: ~34kDa.



Western blot of mouse lung (A), mouse testis (B), rat lung (C) and rat testis (D) lysates with NEK7 antibody at 0.3ug/ml. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: NIMA (never in mitosis gene a)-related kinase 7

Application Notes

Optimal dilution of the NEK7 antibody should be determined by the researcher.

Immunogen

Amino acids CINPDPEKRPDVT were used as the immunogen for this NEK7 antibody.

Storage

Aliquot and store the NEK7 antibody at -20oC.