


NEK7 Antibody (R34933)

Catalog No.	Formulation	Size
R34933-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	140609
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This NEK7 antibody is available for research use only.



Western blot testing of HeLa lysate with NEK7 antibody at 0.05ug/ml. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: NIMA-related kinase 7

Application Notes

Optimal dilution of the NEK7 antibody should be determined by the researcher.

Immunogen

Amino acids NTLANFRIEKK were used as the immunogen for this NEK7 antibody.

Storage

Aliquot and store the NEK7 antibody at -20oC.