


Neuroblastoma Antibody (R34513)

Catalog No.	Formulation	Size
R34513-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4681
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Neuroblastoma antibody is available for research use only.



Western blot testing of KELLY cell lysate with Neuroblastoma antibody at 0.1ug/ml.
Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: Neuroblastoma, suppression of tumorigenicity 1; NBL1

Application Notes

Optimal dilution of the Neuroblastoma antibody should be determined by the researcher.

1. This Neuroblastoma antibody is specific for isoform 1.

Immunogen

Amino acids KLGKTHGHRAL E were used as the immunogen for this Neuroblastoma antibody.

Storage

Aliquot and store the Neuroblastoma antibody at -20oC.