


Neuroblastoma Antibody (R34513)

Catalog No.	Formulation	Size
R34513-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4681
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Neuroblastoma antibody is available for research use only.



Description

Additional name(s) for this target protein: Neuroblastoma, suppression of tumorigenicity 1; NBL1

Application Notes

Optimal dilution of the Neuroblastoma antibody should be determined by the researcher.

1. This Neuroblastoma antibody is specific for isoform 1.

Immunogen

Amino acids KLGKTHGHRALE were used as the immunogen for this Neuroblastoma antibody.

Storage

Aliquot and store the Neuroblastoma antibody at -20oC.