


NeuroD4 Antibody (R34228)

Catalog No.	Formulation	Size
R34228-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8193
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This NeuroD4 antibody is available for research use only.



Western blot testing of human cerebellum lysate with NeuroD4 antibody at 1ug/ml.
Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Zinc and double PHD fingers family 1; DPF1

Application Notes

Optimal dilution of the NeuroD4 antibody should be determined by the researcher.

Immunogen

Amino acids HLKEKASAYITLT were used as the immunogen for this NeuroD4 antibody.

Storage

Aliquot and store the NeuroD4 antibody at -20oC.