


Neurofilament-66 Antibody / Alpha Internexin (RQ6489)


Catalog No.	Formulation	Size
RQ6489	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q16352
Localization	Cytoplasmic
Applications	Western Blot : 1-2ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Immunofluorescence : 5ug/ml Immunofluorescence : 5ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This Neurofilament-66 antibody is available for research use only.




IHC staining of FFPE mouse brain tissue with Neurofilament-66 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE mouse brain tissue with Neurofilament-66 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE rat brain tissue with Neurofilament-66 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




Immunofluorescent staining of FFPE rat brain tissue with Neurofilament-66 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Immunofluorescent staining of FFPE human SiHa cells with Neurofilament-66 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of 1) rat brain and 2) mouse brain tissue lysate with Neurofilament-66 antibody. Predicted molecular weight ~55 kDa but routinely observed at ~66 kDa.



Flow cytometry testing of human A549 cells with Neurofilament-66 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue=Neurofilament-66 antibody.

Description

Alpha-Internexin (INA; also NF-66) is a 66 kDa member of the intermediate filament (IF) protein family. The protein was originally purified from rat optic nerve and spinal cord. And the protein copurifies with other neurofilament subunits, as it was originally discovered, however in some mature neurons it can be the only neurofilament expressed. The protein is present in developing neuroblasts and in the central nervous system of adults. Meanwhile, the protein is a major component of the intermediate filament network in small interneurons and cerebellar granule cells, where it is present in the parallel fibers.

Application Notes

Optimal dilution of the Neurofilament-66 antibody should be determined by the researcher.

Immunogen

An E. coli-derived human protein (amino acids L62-D169) was used as the immunogen for the Neurofilament-66 antibody.

Storage

After reconstitution, the Neurofilament-66 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.