


# Neurogranin Antibody (Biotin Conjugate) (R34473BTN)

Catalog No.	Formulation	Size
R34473BTN	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q92686
Gene ID	4900
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Neurogranin antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: NRGN, formerly designated p17

### **Application Notes**

Optimal dilution of the Neurogranin antibody should be determined by the researcher.

#### **Immunogen**

Amino acids KIKSGERGRKG were used as the immunogen for this Neurogranin antibody.

#### **Storage**

Aliquot and store the Neurogranin antibody at -20oC.

References (1)