


Neurologin 2 Antibody (R34877)

Catalog No.	Formulation	Size
R34877-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57555
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Neurologin 2 antibody is available for research use only.



Western blot testing of human pancreas lysate with Neurologin 2 antibody at 0.1ug/ml.
Predicted molecular weight: ~91kDa.

Description

Additional name(s) for this target protein: NLGN2

Application Notes

Optimal dilution of the Neurologin 2 antibody should be determined by the researcher.

Immunogen

Amino acids NPPD⁺TD⁺IRDPGKK were used as the immunogen for this Neuroligin 2 antibody.

Storage

Aliquot and store the Neuroligin 2 antibody at -20oC.