


Neurotrophin 4 Antibody (R34057)

Catalog No.	Formulation	Size
R34057-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4909
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Neurotrophin 4 antibody is available for research use only.



Western blot testing of human cerebellum lysate with Neurotrophin 4 antibody at 0.1ug/ml. Expected molecular weight: 22-45 depending on level of glycosylation.

Description

Additional name(s) for this target protein: NT5 antibody, NTF4, neurotrophic factor 4, neurotrophic factor 5; NTF5

Application Notes

Optimal dilution of the Neurotrophin 4 antibody should be determined by the researcher.

Immunogen

Amino acids FETRCKADNAEEG were used as the immunogen for this Neurotrophin 4 antibody.

Storage

Aliquot and store the Neurotrophin 4 antibody at -20oC.

References (1)