


Neurturin Antibody (R33378)

Catalog No.	Formulation	Size
R33378-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4902
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Neurturin antibody is available for research use only.



Western blot testing of human heart lysate with Neurturin antibody at 0.1ug/ml. Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: NTN; NRTN

Application Notes

Optimal dilution of the Neurturin antibody should be determined by the researcher.

Immunogen

Amino acids DAHSRYHTVHELARE were used as the immunogen for this Neurturin antibody.

Storage

Aliquot and store the Neurturin antibody at -20oC.