


## Neurturin Antibody (R33667)

Catalog No.	Formulation	Size
R33667-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	18188
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Neurturin antibody is available for research use only.



Western blot testing of mouse brain lysate with Neurturin antibody at 0.1ug/ml. The expected ~23kDa band and the additional ~50kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: NRTN; NTN

## Application Notes

Optimal dilution of the Neurturin antibody should be determined by the researcher.

## Immunogen

Amino acids DVHSRYHTLQELSARE were used as the immunogen for this Neurturin antibody.

## **Storage**

Aliquot and store the Neurturin antibody at -20oC.