


# **Neurturin Antibody (R33667)**

Catalog No.	Formulation	Size
R33667-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	18188
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Neurturin antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: NRTN; NTN

### **Application Notes**

Optimal dilution of the Neurturin antibody should be determined by the researcher.

### **Immunogen**

Amino acids DVHSRYHTLQELSARE were used as the immunogen for this Neurturin antibody.

# **Storage** Aliquot and store the Neurturin antibody at -20oC.