


NFATC3 Antibody / NFAT4 (R34255)

Catalog No.	Formulation	Size
R34255-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4775
Applications	Western Blot : 0.3-1.0ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This NFATC3 antibody is available for research use only.



Western blot testing of Jurkat nuclear lysate with NFATC3 antibody at 0.03ug/ml. Predicted molecular weight: ~116 kDa but can be observed at 140~150 kDa due to phosphorylation.

Description

Additional name(s) for this target protein: Nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3, NFAT4

Application Notes

Optimal dilution of the NFATC3 antibody should be determined by the researcher.

Immunogen

Amino acids HQELDAHEDDLQIN were used as the immunogen for this NFATC3 antibody.

Storage

Aliquot and store the NFATC3 antibody at -20oC.