


NFIL3 Antibody (R35669)

Catalog No.	Formulation	Size
R35669-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4783
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This NFIL3 antibody is available for research use only.



Western blot testing of human Daudi lysate with NFIL3 antibody at 2ug/ml. Predicted molecular weight: ~52 kDa.

Description

Additional name(s) for this target protein: Nuclear factor interleukin 3 regulated

Application Notes

Optimal dilution of the NFIL3 antibody should be determined by the researcher.

Immunogen

Amino acids RLIATQPISASDSG were used as the immunogen for this NFIL3 antibody.

Storage

Aliquot and store the NFIL3 antibody at -20oC.