


## NGAL Antibody (R34280)

Catalog No.	Formulation	Size
R34280-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16819
Localization	Cytoplasmic
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This NGAL antibody is available for research use only.



Western blot testing of mouse kidney lysate with NGAL antibody at 2ug/ml. Predicted molecular weight: 22-25kDa.

## Description

Additional name(s) for this target protein: Lipocalin 2; Lcn2

## Application Notes

Optimal dilution of the NGAL antibody should be determined by the researcher.

## Immunogen

Amino acids HRYPQVQSYNVQ were used as the immunogen for this NGAL antibody.

## Storage

Aliquot and store the NGAL antibody at -20oC.