


NGF Antibody / Beta Nerve Growth Factor (R30235)


Catalog No.	Formulation	Size
R30235	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P01138
Applications	Western Blot : 0.5-1ug/ml
Limitations	This NGF antibody is available for research use only.



Western blot testing of NGF antibody and rat brain tissue lysate



Western blot testing of NGF antibody and Lane 1: recombinant human protein 10ng; 2: 5ng

Description

Nerve growth factor is a polypeptide involved in the regulation of growth and differentiation of sympathetic and certain sensory neurons. The nucleotide sequence of human and mouse beta-NGF are very similar. The beta-subunits of nerve growth factor(NGFB) have been assigned to mouse chromosome 3 and human chromosome 1p22. The human gene for the beta subunit of nerve growth factor is located on the proximal short arm of chromosome 1. A mutation in the nerve growth factor beta gene(NGFB) causes loss of pain perception.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the NGF antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the N-terminus of human NGF (EVNINNSVFKQY) was used as the immunogen for this NGF antibody.

Storage

After reconstitution, the NGF antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.