


NGFR Antibody (R33143)

Catalog No.	Formulation	Size
R33143-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4804
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This NGFR antibody is available for research use only.



Western blot testing of human heart lysate with NGFR antibody at 0.3ug/ml. Predicted molecular weight ~45/75kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Nerve growth factor receptor, CD271, p75NTR

Application Notes

Optimal dilution of the NGFR antibody should be determined by the researcher.

Immunogen

Amino acids AKREEVEKLLNG were used as the immunogen for this NGFR antibody.

Storage

Aliquot and store the NGFR antibody at -20oC.