


NgR3 Antibody / Nogo-66 receptor-related protein 3 / RTN4RL1 [clone 29R59] (RQ8765)

Catalog No.	Formulation	Size
RQ8765	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	29R59
Purity	Affinity chromatography
UniProt	Q86UN2
Applications	Western Blot : 1:500-1:2000
Limitations	This NgR3 antibody is available for research use only.



Western blot testing of human K562 cell lysate with NgR3 antibody. Predicted molecular weight ~49 kDa.

Description

Like other family members, NgR3/RTN4RL1 plays a role in restricting the number dendritic spines and the number of synapses that are formed during brain development. [UniProt]

Application Notes

Optimal dilution of the NgR3 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human NgR3 protein was used as the immunogen for the NgR3 Antibody.

Storage

Store the NgR3 antibody at -20oC.