


NKX1.1 Antibody (F41779)


Catalog No.	Formulation	Size
F41779-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41779-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q15270
Applications	Western Blot : 1:1000 Immunofluorescence : 1:10-1:50
Limitations	This NKX1.1 antibody is available for research use only.



NKX1.1 antibody western blot analysis in SK-N-MC lysate.



Fluorescent confocal image of HeLa cell stained with NKX1.1 antibody. For secondary Ab, Alexa Fluor 488 conjugated donkey anti-rabbit Ab (green) was used (1:400, 50 min at 37°C). NKX1.1 immunoreactivity is localized to the nucleus.

Description

NKX1.1 belongs to the NK-1 homeobox family. It is a DNA dependent regulator of transcription.

Application Notes

Titration of the NKX1.1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-1 from the human protein was used as the immunogen for this NKX1.1 antibody.

Storage

Aliquot the NKX1.1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.