

NKX2-1 Antibody / NK2 homeobox 1 [clone MSVA-312R] (V6121)

Catalog No.	Formulation	Size
V6121-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6121-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	MSVA-312R
UniProt	P43699
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This NKX2-1/NK2 homeobox 1 antibody is available for research use only.

Application Notes

1. Optimal dilution of the NKX2-1/NK2 homeobox 1 antibody should be determined by the researcher.
2. This NKX2-1/NK2 homeobox 1 antibody is recombinantly produced by expression in CHO cells.

Immunogen

Recombinant full-length human NKX2-1 (TTF-1) protein was used as the immunogen for the NKX2-1/NK2 homeobox 1 antibody.

Storage

NKX2-1/NK2 homeobox 1 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

