


NKX2.4 Antibody (F43599)

Catalog No.	Formulation	Size
F43599-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43599-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9H2Z4
Applications	Western Blot : 1:1000
Limitations	This NKX2.4 antibody is available for research use only.



Description

NKX2-4 is a probable transcription factor.

Application Notes

Titration of the NKX2.4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 2-4 from the human protein was used as the immunogen for this NKX2.4 antibody.

Storage

Aliquot the NKX2.4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.