


NKX2.6 Antibody (F46308)

Catalog No.	Formulation	Size
F46308-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F46308-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	A6NCS4
Applications	Western Blot : 1:1000
Limitations	This NKX2.6 antibody is available for research use only.



Description

In conjunction with NKX2-5, may play a role in both pharyngeal and cardiac embryonic development (By similarity).

Application Notes

Titration of the NKX2.6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 2-6 from the human protein was used as the immunogen for this NKX2.6 antibody.

Storage

Aliquot the NKX2.6 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.