


## NKX3.2 Antibody (F45860)

Catalog No.	Formulation	Size
F45860-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F45860-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Mouse, Xenopus
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P78367
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This NKX3.2 antibody is available for research use only.



NKX3.2 antibody western blot analysis in U266 lysate.

## Description

This gene encodes a member of the NK family of homeobox-containing proteins. The encoded protein may play a role in skeletal development.

## Application Notes

Titration of the NKX3.2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

A portion of amino acids 3-2 from the human protein was used as the immunogen for this NKX3.2 antibody.

## **Storage**

Aliquot the NKX3.2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.