


NLGN4X Antibody (R34949)

Catalog No.	Formulation	Size
R34949-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57502
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NLGN4X antibody is available for research use only.



Description

Additional name(s) for this target protein: Neuroligin 4 X-linked

Application Notes

Optimal dilution of the NLGN4X antibody should be determined by the researcher.

Immunogen

Amino acids NYGKIRGLRTPLPNE were used as the immunogen for this NLGN4X antibody. **Storage** Aliquot and store the NLGN4X antibody at -20oC.