


NLGN4X Antibody (R34949)

Catalog No.	Formulation	Size
R34949-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57502
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NLGN4X antibody is available for research use only.



Western blot testing of amygdala lysate with NLGN4X antibody at 1ug/ml. Predicted molecular weight: ~92kDa.

Description

Additional name(s) for this target protein: Neuroligin 4 X-linked

Application Notes

Optimal dilution of the NLGN4X antibody should be determined by the researcher.

Immunogen

Amino acids NYGKIRGLRTPNE were used as the immunogen for this NLGN4X antibody.

Storage

Aliquot and store the NLGN4X antibody at -20oC.