


Nlrp12 Antibody (R34567)

Catalog No.	Formulation	Size
R34567-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	378425
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Nlrp12 antibody is available for research use only.



Western blot testing of mouse spleen lysate with Nlrp12 antibody at 1ug/ml. Predicted molecular weight: ~119kDa, routinely observed at 100~119kDa.

Description

Additional name(s) for this target protein: NACHT, LRR and PYD containing protein 12, Nalp12

Application Notes

Optimal dilution of the Nlrp12 antibody should be determined by the researcher.

Immunogen

Amino acids KPGTPTFKVPANQR were used as the immunogen for this Nlrp12 antibody.

Storage

Aliquot and store the Nlrp12 antibody at -20oC.