


Nmnat3 Antibody (R34352)

Catalog No.	Formulation	Size
R34352-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	74080
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Nmnat3 antibody is available for research use only.



Western blot testing of mouse heart lysate with Nmnat3 antibody at 0.5ug/ml. Predicted/observed molecular weight: ~28kDa. Both observed bands are blocked by addition of immunizing peptide.

Description

Additional name(s) for this target protein: Nicotinamide nucleotide adenylyltransferase 3

Application Notes

Optimal dilution of the Nmnat3 antibody should be determined by the researcher.

Immunogen

Amino acids KALGQGQSVKYLLPE were used as the immunogen for this Nmnat3 antibody.

Storage

Aliquot and store the Nmnat3 antibody at -20oC.