


NNMT Antibody (R34976)

Catalog No.	Formulation	Size
R34976-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4837
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This NNMT antibody is available for research use only.



Western blot testing of HeLa lysate with NNMT antibody at 0.3ug/ml. Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: Nicotinamide N-methyltransferase

Application Notes

Optimal dilution of the NNMT antibody should be determined by the researcher.

Immunogen

Amino acids PDLPTYCRALRN were used as the immunogen for this NNMT antibody.

Storage

Aliquot and store the NNMT antibody at -20oC.