


Nogo 66 receptor Antibody (R34733)

Catalog No.	Formulation	Size
R34733-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	65078
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Nogo 66 receptor antibody is available for research use only.



Western blot testing of human brain lysate with Nogo 66 receptor antibody at 0.3ug/ml with [B] and without [A] blocking/immunizing peptide

Description

Additional name(s) for this target protein: Reticulon 4 receptor; RTN4R

Application Notes

Optimal dilution of the Nogo 66 receptor antibody should be determined by the researcher.

Immunogen

Amino acids LQALPDDTFRDLGN were used as the immunogen for this Nogo 66 receptor antibody.

Storage

Aliquot and store the Nogo 66 receptor antibody at -20oC.