


NOXO1 Antibody (R35614)

Catalog No.	Formulation	Size
R35614-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	124056
Applications	Immunohistochemistry (FFPE) : 10-20ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This NOXO1 antibody is available for research use only.



IHC testing of FFPE human kidney with NOXO1 antibody at 10ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: NADPH oxidase organizer 1

Application Notes

Optimal dilution of the NOXO1 antibody should be determined by the researcher.

Immunogen

Amino acids RGCVDSVPHPTTEQ were used as the immunogen for this NOXO1 antibody.

Storage

Aliquot and store the NOXO1 antibody at -20oC.