


## NPC2 Antibody / Niemann Pick C2 (R35186)

Catalog No.	Formulation	Size
R35186-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10577
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This NPC2 antibody is available for research use only.



Western blot testing of human testis lysate with NPC2 antibody at 0.5ug/ml. Predicted molecular weight: ~17 kDa, can be observed as a ~21/23 kDa doublet in some samples. (Ref 1).

### Description

Additional name(s) for this target protein: Niemann-Pick disease, type C2

### Application Notes

Optimal dilution of the NPC2 antibody should be determined by the researcher.

### Immunogen

Amino acids QAEPVQFKDCGSVD were used as the immunogen for this NPC2 antibody.

## Storage

Aliquot and store the NPC2 antibody at -20oC.

## References (1)