


Nprl3 Antibody (R35443)

Catalog No.	Formulation	Size
R35443-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8131
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Nprl3 antibody is available for research use only.



Western blot testing of human heart lysate with Nprl3 antibody at 0.3ug/ml. Predicted molecular weight: ~63kDa, observed here at ~55kDa.

Description

Additional name(s) for this target protein: Nitrogen permease regulator-like 3

Application Notes

Optimal dilution of the Nprl3 antibody should be determined by the researcher.

Immunogen

Amino acids QRSQEHPASQTNKP were used as the immunogen for this Nprl3 antibody.

Storage

Aliquot and store the Nprl3 antibody at -20oC.