


# NPTX1 Antibody (R33983)

| Catalog No.  | Formulation  | Size   |
|--------------|--|--------|
| R33983-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$ | 100 ug |

#### **Bulk quote request**

| Availability       | 1-3 business days                                       |
|--------------------|---|
| Species Reactivity | Human   |
| Format             | Antigen affinity purified                               |
| Clonality          | Polyclonal (goat origin)                                |
| Isotype            | Goat Ig   |
| Purity             | Antigen affinity  |
| Gene ID            | 4884  |
| Applications       | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000  |
| Limitations        | This NPTX1 antibody is available for research use only. |



## **Description**

Additional name(s) for this target protein: Neuronal pentraxin 1

## **Application Notes**

Optimal dilution of the NPTX1 antibody should be determined by the researcher.

### **Immunogen**

Amino acids ERVKIETALTSLHQR were used as the immunogen for this NPTX1 antibody.

| Storage Aliquot and store the NPTX1 antibody at -20oC. |  |  |  |  |
|--|--|--|--|--|
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |
|  |  |  |  |  |